




LSVT Companion[®] Home Edition

LSVT Global, Inc. created the Lee Silverman Voice Treatment (LSVT LOUD[®]) to help individuals with Parkinson disease and other neurological conditions regain and maintain their speaking voice. LSVT LOUD has helped tens of thousands of individuals worldwide with effective, long-lasting speech therapy.

LSVT Global developed the LSVT Companion System, Home Edition to help individuals practice their LOUD voice exercises at home during and after treatment with an LSVT LOUD Certified Clinician. The LSVT Companion, Home Edition is available to individuals working with an LSVT LOUD Certified Clinician and is the only solution of its type cleared by the FDA for use as a medical device.

Your Personal Treatment Plan

The LSVT Companion, Home Edition supports your personal treatment plan created by your speech-language clinician. Using their LSVT Companion, Clinician Edition, your speech-language clinician personalizes your speech exercises and sets goals for you by defining target performance levels and number of repetitions for each LSVT LOUD voice exercise. As you progress, your speech-language clinician can easily update your goals.



LSVT Companion displays client performance during voice and speech exercises.

Convenient and Effective Practice

Practicing your LOUD exercises at home will augment in-person sessions with your speech-language clinician during treatment and help maintain your LOUD voice after you complete the in-person sessions. The voice and speech exercises on the LSVT Companion will be familiar to you from having done the same LOUD exercises in your speech-language clinician's office using their LSVT Companion.

To help you receive the most benefit from your home voice exercises, as you perform each voice exercise the LSVT Companion displays key measures including loudness (sound pressure level), frequency, and duration. In addition, the LSVT Companion has a "virtual coach" that assesses how your performance compares with goals set by your speech-language clinician and provides audio feedback through speakers or headphones attached to your computer.

Supporting Your Clinician Partnership

Just like using the LSVT Companion in your speech-language clinician's office, the LSVT Companion, Home Edition captures your vocal performance data for each of your voice and speech exercises and stores those results on your computer. You can easily provide those data to your speech-language clinician so they can evaluate your progress and compare your home practice with in-person treatment sessions.

The LSVT Companion, Home Edition includes a microphone calibrated using LSVT Global patent-pending process to collect accurate voice exercise data that is consistent with the data collected in your speech-language clinician's office.





LSVT Companion[®] Home Edition

LSVT Companion, Home Edition Licensing

The LSVT Companion, Home Edition is licensed for use by a single user. Included with the LSVT Companion, Home Edition is one calibrated microphone and the software program.

The LSVT Companion, Home Edition includes software and a calibrated microphone for collecting accurate data. The LSVT Companion, Home Edition is licensed for use by a single user. A client must be working with and LSVT LOUD Certified Clinician who has purchased the LSVT Companion, Clinician Edition to use this product.

The LSVT Companion is cleared for use as a medical device by the FDA and has CE registration. At this time the LSVT Companion, Home Edition can be sold only in the U.S., Canada, UK, Ireland, Cypress, Estonia, Finland, France, Greece, Iceland, Luxembourg, Malta, Poland, New Zealand, and Australia.

Minimum System Requirements

- Microsoft[®] Compatible 1GHz Pentium Class Computer
- 512 Mb RAM
- Windows[®] 2000 SP4, XP[®] SP2, XP Professional[®] SP2, VISTA[®], Windows[®] 7, or Windows[®] 8 operating systems
- Microsoft Excel (for viewing data)
- Headphones or computer speakers

The LSVT Companion, Home Edition makes it easy for individuals to practice their LSVT LOUD exercises at home, helping them regain and maintain their LOUD voice.
Lorraine Ramig, Ph.D., CCC-SLP

If you have any questions about the LSVT Companion, we look forward to hearing from you!

LSVT Global, Inc.

**3323 N Campbell Avenue, Suite 5
Tucson, AZ 85719 USA**

Toll Free: 1-888-438-5788 (in the USA)

Phone: 1-520-867-8838

Email: info@LSVTGlobal.com

Website: www.LSVTGlobal.com

Fax: 520-867-8839





LSVT Companion[®] Home Edition

ITEMS

	Price	Quantity	Total
LSVT Companion, Home Edition	\$299.00 each	_____	_____
Sales Tax - <i>Arizona orders only, 8.1%</i>	\$24.22 each	_____	_____
Shipping and handling - <i>Continental U.S.*</i>	\$13.45 each	_____	_____
<i>*USPS[®] international flat rates may vary. Contact us for further assistance.</i>			
		TOTAL:	_____

Name of the LSVT LOUD Certified Clinician with whom the client is working:

METHOD OF PAYMENT *Choose a payment option:*

Check

Master Card

VISA

AMEX

Card #: _____

Expiration Date: _____

Cardholder Name: _____

Security Code: _____

BILLING INFORMATION (required)

If paying by credit card, please use billing address associated with the card.

First Name: _____ Middle Initial: _____

Last Name: _____

Affiliation: _____

Address: _____

Apt/Suite/Other: _____

City: _____ State: _____ Zip Code: _____

Phone #: _____ Email: _____

SHIPPING INFORMATION

Check if Shipping and Billing information is the same.

First Name: _____ Middle Initial: _____

Last Name: _____

Affiliation: _____

Address: _____

Apt/Suite/Other: _____

City: _____ State: _____ Zip Code: _____

Phone #: _____ Email: _____

Please complete this form and submit your order using one of the following methods:

- MAIL: **LSVT Global, 3323 N Campbell Ave, Suite 5, Tucson, AZ 85719**
- FAX: **(520) 867-8839**
- EMAIL: **info@LSVTGlobal.com**
- ONLINE: **www.LSVTGlobal.com**

